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Instruction For Use

Urine Bag(Sterile)

Type: Leg Bag, Drainage Bag (Type A, Type B)

DESCRIPTION

The Urine Bag is used for collection of urine acting as a reservoir attached to the end of a urinary catheter. The product is used in hospitals or in home health care with patients using Bladder- or Male External Catheters.

FEATURES

- Can be used in hospitals and in home situations where an individual may be suffering from chronic illness or conditions that make regular trips to the bathroom unfeasible.
- Provides an efficient solution for urinary incontinence. A leg bag preserves freedom of mobility and activity as the bag is attached to the leg of the patient.

INTENDED USE

The Urine Bag is used for the collection and temporary storage of urine. It is intended to be used with an indwelling catheter in incontinent patients, cannot urinate normally, or are requiring continuous bladder drainage.

CLINICAL BENEFITS

By connecting to the end of the urinary catheter or male external catheter, the Urine Bag can be used for the collection and temporary storage of urine. It is used for patients suffering from chronic illness or conditions that make regular trips to the bathroom unfeasible. The clinical benefit associated with the device is define as direct, it is intended to maintain a person's mobility, dignity, and comfort.

PATIENT TARGET GROUP

Neonate, Infant, Child, Adult.

INTENDED USER

The Urine Bag can be operated by trained professionals or lay persons who are trained on the drainage procedure. If operated by lay persons, any question during operation or use, consult trained professionals.

INDICATIONS

The Urine Bag is suitable for connecting to a Foley (bladder) catheter for measuring patient urine and urine collection.

CONTRAINDICATIONS

There are no known contraindications.

ADVERSE REACTIONS

No adverse reaction has been reported.

DESCRIPTION FOR USE

Leg Bag

- Open the package and take out the Leg Bag.
- Fix the Leg Bag on the patient's leg.
- Connect the connector to the drainage funnel of the Catheter.
- When the urine reaches the upper limit of the leg bag, drain the urine from bag through the bottom outlet without disconnecting the hag from the catheter, if necessary, exchange the leg bag with a new one.
- Open the bottom outlet to drain the urine into the toilet or a Suitable hygienic disposal. If required record the amount of collected urine.

Drainage Bag (Type A, Type B)

- Open the package and take out the Urine Bag from the package.
- Hang the Bag near the foot of the bed below bladder level and position the connecting tube to allow an easy flow. Connect the connector to the drainage funnel of the Catheter.
- When the urine reaches the upper limit of the Bag, drain the urine from Drainage the Bag through the bottom outlet (if available) without disconnecting the bag from the catheter. If necessary, exchange the Bag with a new one.
- Open the bottom outlet (if available) to drain the urine into the toilet or a suitable hygienic disposal. If required record the amount of collected urine.

Exchange of new Leg Bag /Drainage Bag (Type A, Type B)

- Disconnect the bag from the catheter.
- Close the tube clip (if available).
- Clean the drainage funnel of the Foley Catheter.
- Connect other Bag to the Catheter following the instructions above.

Janmed+

Instruction For Use

PRECAUTIONS/WARNINGS

PRECAUTIONS

- The Urine Bag should be placed below bladder level. If used in hospital, follow the hospital procedures to use the product. If used in home, consult health care professionals for professional directions.
- Do not kink or squeeze the tubing, as this will affect the drainage of urine.
- Do not place the bag on the bed. The bag must be positioned in a hanging position, otherwise there is a risk of the bag leaking.
- Duration of use - Max 7 days.

WARNINGS

- Read all warnings and instructions before use. Improper use can result in serious or fatal illness or injury.
- Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Do not use if it's damaged or irregularly shaped.
- Do not use after the expiry date.
- For single patient use only.
- Dispose of product and packaging in accordance with hospital administrative and/or local government policy.

STORAGE CONDITIONS

- Store in a cool and dry place, keep away from sunlight.
- Protect product from moisture and excessive heat.
- Avoid prolonged exposure to ultraviolet, sunlight and fluorescent light.
- Store in manner preventing crushing.

The Key Performance Characteristics of Urine Bag

- Freedom from leakage under impact: water is added to the Urine Bag, which is then dropped on to a smooth, rigid surface, the bag must not leak.

MEANING OF SYMBOLS ON PACKAGE

	Date of manufacture		Use-by date		Batch code
	Do not re-use		Consult instructions for use		Do not use if package is damaged
	Sterilized using ethylene oxide		Do not re-sterilize		Single sterile barrier system
	Medical Device		Manufacturer		Authorized representative in the European Community
	CE Marked Product		Keep dry		Keep away from sunlight
	Importer		This way up		Fragile, handle with care
	Catalogue number		Temperature limits +15°C -49°C		UKCA Marked Product
	UK Responsible Person		Authorized representative in Switzerland		Unique device identifier
	Does not contain or presence of phthalate				

Janmed GmbH, Otto-Hahn-Straße 18, 71069 Sindelfingen, Germany
Tel.: +49(0)7031-76234-0, Fax: +49(0)7031-76234-19, info@janmed.de, www.janmed.de

WELL LEAD MEDICAL CO., LTD.
Address: C-4 Jinhui Industrial Estate, Huakong, 511434 Panyu, Guangzhou, PEOPLE'S REPUBLIC OF CHINA
E-mail: info@welllead.com.cn www.welllead.com.cn
Tel: +86 20 84758878

Shanghai International Holding Corp. GmbH (Europe)
Address: Effelstrasse 80, 20537 Hamburg, GERMANY

SUNGO Certification Company Limited
Address: 3rd Floor, 70 Gracechurch Street, London, EC3V 9HR England, United Kingdom

QS Engineering AG
Address: Edenstrasse 31, CH - 4106 Therwil, Switzerland

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Well Lead Medical CO.,LTD